



Quick Start Guide: How to Install and Start 3D Slicer

Sonia Pujol, Ph.D.

Department of Radiology Brigham and Women's Hospital Harvard Medical School


spujol@bwh.harvard.edu



Step 1:

- Go the 3D Slicer website: <u>www.slicer.org</u>
- Click on Download





Step 2:


ullet

Choose the stable
release version
appropriate for your
platform and start
downloading the
software




Step 3:


 Once the download is complete, doubleclick on the Slicer package to start the installing application



Step 4:

 After installing the application, doubleclick on the Slicer icon to start running the software






Note for Mac Users:

- A pop-up window might indicate the Slicer application cannot be installed
- Go to Security Preferences and select Security & Privacy

	Contractor Contractor Contractor	Change Password	
Require password	1 hour 📀	after sleep or screen say	er begins
Show a message	when the screen is I	ocked Set Lock Messa	ge
Allow apps downloaded f	rom:		
Allow apps downloaded f	rom:		
App Store	entified developers		Open Anywa

Note for Mac Users:

- Select Allow apps downloaded from the App Store and identified developers
- Click on **Open Anyway**



3D Slicer release version 4

spujol@bwh.harvard.edu